



PERSONAL PROTECTIVE EQUIPMENT (PPE)

**COMPLIES WITH ANSI Z87.1 CSA Z94.3 EN166
AND AS/NZ 1337 STANDARDS**

6100 Hollister Avenue
Building 3, Unit 3A/3B
Santa Barbara, CA 93117

P: 805.967.1153
F: 805.967.1295

tecfen@tecfen.com
www.tecfenmedical.com
BR10_Rev1



Face Shield

- Anti-Fog
- Latex Free
- PET material 0.188mm thick
- Adjustable strap for universal fit
- 380mm width x 230mm height for maximum coverage



Face Shield S

- Anti-Fog
- Latex Free
- PET material 0.188mm thick
- Hooks on end of frame allow for addition of rubber/elastic band for increased hold
- 230 width x 200mm height for excellent coverage




PRODUCT NO.	DESCRIPTION	QTY
FS8001	159mm W x 30mm D x 230mm H	200/Case

PRODUCT NO.	DESCRIPTION	QTY
FS8002	159mm W x 179mm D x 200mm H	200/Case

Face Shield XS

- Anti-Fog
- Latex Free
- PET material 0.188mm thick
- Hooks on end of frame allow for addition of rubber/elastic band for increased hold
- 230mm width x 69mm height for streamlined coverage



PRODUCT NO.	DESCRIPTION	QTY
FS8003	157mm W x 160mm D x 69mm H	240/Case

Goggles

- Anti-Fog
- Latex Free
- Clear PVC material forms seal around face
- Fits over most prescription eyewear
- 155mm width x 75mm height



PRODUCT NO.	DESCRIPTION	QTY
G200S	Goggles with Indirect Ventilation 155mm W x 75mm H x 71mm D	120/Case
G200CS	Goggles without Ventilation 155mm W x 75mm H x 71mm D	120/Case

Protective glasses

- Eyewear lens and temple is made of Polycarbonate PC
- Latex Free
- The frame contains no metallic parts
- Fits over most prescription eyewear



PRODUCT NO.	DESCRIPTION	QTY PER CASE
Y20C	Protective Glasses (149mm W x 57mm H)	120/Case