



LASER PROBES


**COMPATIBLE WITH ALCON, IRIDEX, COHERENT,
ZEISS, ELLEX, QUANTEL AND DORC LASERS**

6100 Hollister Avenue
Building 3, Unit 3A/3B
Santa Barbara, CA 93117

P: 805.967.1153
F: 805.967.1295

tecfen@tecfen.com
www.tecfenmedical.com
BR8_Rev1





While there are many similarities between laser probes these days, Our line focuses on several key performance areas to deliver a significant improvement to your endo-photocoagulation experience.

- Smooth extension and retraction of curve.
- Uniform spot size when straight or curved.
- Consistant laser uptake.
- Direct access to posterior pole.
- Fixed flexible curve allows easy access through rigid cannula.
- Smooth transition at taper.

Straight Laser Probe:

Direct Access to posterior pole, uniform spot size



PRODUCT NO.	DESCRIPTION	QTY PER CASE
QLPS-23	Straight Laser Probe, 23g	5/box
QLPS-25	Straight Laser Probe, 25g	5/box
QLPS-26	Straight Laser Probe, 27g	5/box

Curved Laser Probe:

Curve angle of 34° for easier access to peripheral retina



PRODUCT NO.	DESCRIPTION	QTY PER CASE
QLPC-23	Curved Laser Probe, 23g	5/box
QLPC-25	Curved Laser Probe, 25g	5/box
QLPC-26	Curved Laser Probe, 27g	5/box

Flexible Laser Probe:

Fixed flexible curve allows easy access through rigid cannula



PRODUCT NO.	DESCRIPTION	QTY PER CASE
QLPF-23	Flexible Laser Probe, 23g	5/box
QLPF-25	Flexible Laser Probe, 25g	5/box

Adjustable Laser Probe:

Smooth extension and retraction of curve



PRODUCT NO.	DESCRIPTION	QTY PER CASE
QLPA-23	Adjustable Laser Probe, 23g	5/box
QLPA-25	Adjustable Laser Probe, 25g	5/box