



# MEMBRANE POLISHERS

USED FOR NEUROEPITHELIAL  
MEMBRANE POLISHING DURING  
RETINAL SURGERY

- Black silicone tip maximizes the contrast against the retina.

- Diamond dust ensures the least amount of shed particles in the eye.



**Membrane Polisher:** .....

PRODUCT NO.	DESCRIPTION	QTY PER CASE
QTP1267-23	Membrane Polisher, 23g	5/box
QTP1267-25	Membrane Polisher, 25g	5/box



**Retractable Membrane Polisher:** .....

PRODUCT NO.	DESCRIPTION	QTY PER CASE
QTPR1267-23	Retractable Membrane Polisher, 23g	5/box
QTPR1267-25	Retractable Membrane Polisher, 25g	5/box